

Avaco

Vacuum therapy



Features

product code A-UV-AST-VCWH operation mosynchronized increased low independent treatment channels

A-UV-AST-VCWH operation mosynchronized increased low increased low increased low in single-charges self-sealing sel

Vacuum therapy

operation modes: continuous, pulse, variable pulse, synchronized with the electric current, pulse with increased lower vacuum

two independent treatment channels (easy operation in single-channel mode)

self-sealing suction cups

possibility of vacuum adjustment

current parameter control during the treatment procedure

container filling signalization

dedicated mode for emptying the water container

independent sockets for the application of electrotherapy electrodes

Vacuum therapy technical parameters

vacuum wave pulse frequency 6, 12, 20, 40, 60 imp/min
vacuum maks 500 mbar/hPa
vacuum adjustment range 100 - 500 mbar, step 40
mbar

General technical parameters

dimensions 30 x 26 x 12 cm

device weight 3,5 kg

power supply, power consumption 230 V, 50 Hz, 20 W, 40





Avaco page 1 / 1