

MAGNUM L

MAGNETOTHERAPY | 8 PROGRAMS

PHYSIOTHERAPY PRODUCT SUITABLE TO RECOVER THE FUNCTIONALITY COMPROMISED BY TRAUMA.

8 HEALTH AND REHAB

TECHNICAL FEATURES

Display

Visible area size 2,6"

Channels

1

Emitted intensity

120 average Gauss

Adjustable frequency

5-200Hz

Duty Cycle

5-50%

Solenoids

1 flexible size mm 300x100

Power supply

Mains

Size


mm 100x160x35

Weight

gr 280

EQUIPMENT

- 1 Bag
- 1 MAGNUM L magnetotherapy unit
- 1 Flexible solenoid mm 300x100
- 1 Power supply unit
- 1 Operating manual



CE
0476



LAST 10



FAVORITES



STIM LOCK



MAINS SUPPLY

GLOBUS
ITALIAN EXCELLENCE

MAGNUM L

MAGNETOTHERAPY | 8 PROGRAMS

PROGRAMS LIST

HEALTH and REHAB

Small fractures

Big fractures

Cervical Osteoarthritis

Back Pain

Osteoporosis

Arthrosis

Arthritis

Epicondylitis

